



Interfacial Shear Strength (IFSS) Module for LEX820



The Dia-Stron Interfacial Shear Strength module (IFSS) is an interchangeable module for the LEX820 high resolution extensometer used to measure the debonding force of micro-droplets on single filaments and fibres.


General Information

Principal Features

- 50mm linear travel
- Highly accurate speed control
- 2.5N & 20N load cells available
- Standard set of shearing plates

Principal Benefits

- Exceptionally smooth travel
- High positional repeatability
- Highly detailed debonding data



Above: Lyocell fiber with polypropylene droplet
(Courtesy of Hochschule Bremen)

Introduction

Fibre-matrix interfacial properties are critical to achieving satisfactory composite material performance. The IFSS module is an interchangeable accessory designed to measure the debonding force of micro-droplets on single filaments and fibres.

This IFSS measurement is based on the universally recognised micro-bond method to evaluate the interfacial properties of a variety of matrix resins and epoxies on fibers and filaments commonly used in composite materials. The IFSS method can be applied to various fibre and filament types: glass, carbon, ceramic, aramid, basalt or natural fibres.

Delivering Measurement Solutions

e: sales@diastron.com

w: www.diastron.com


UK office 9 Focus Way | Andover | Hampshire | SP10 5NY | UK
t: +44 (0)1264 334700 | f: +44 (0)1264 334686

888 Sussex Boulevard | Broomall | PA19008 | USA
t: (610) 328-9038 | f: (610) 328-7157

US office

The company reserves the right to alter specifications and prices without giving prior notice.

© Dia-Stron Limited



Specifications

LEX820

Extension range	3 - 53mm
Speed range	0.01 to 2.6mm/sec
Force range	0 to 2.5N or 0 to 20N
Force resolution	0.05mN (2.5N) 0.5mN (20N)
Displacement resolution	1µm
Displacement accuracy	50µm
Load cell linearity	±0.1% full scale

Shearing plates

Standard shearing plate slot sizes	<ul style="list-style-type: none"> ● 50µm ● 80 µm ● 100 µm ● 200 µm
------------------------------------	---

Content


LEX820 Instrument
IFSS Module
UV1000 Control unit
PU1100 Pneumatic unit
UvWin software for Windows OS

Requirements

Power Supply	85-265vac 47-63Hz, 100W
Compressed Air: min. 5 Bar	
Computer	<ul style="list-style-type: none"> ● Windows OS: 7, 8, 10 ● 1 x USB port


IFSS Module

The IFSS module uses interchangeable precision laser cut tungsten plates to support the micro-droplet whilst the specimen is withdrawn through. The force being applied by the micro-droplet to the plate is recorded by the load cell until interfacial failure. The sample is secured at the other end using the Dia-Stron one part plastic tab system and held in place using pneumatic sample covers.



Dedicated software – UvWin

UvWin 4 software controls the IFSS system. Method parameters can be easily edited within the software. UvWin enables automatic correction for system compliance.



Debonding data for a polypropylene droplet from a Lyocell fibre

UvWin also offers a number of integrated data processing tools to analyse the data. The raw data can also be exported.

Sample Mounting

Samples are mounted using the Dia-Stron one part plastic tab system. Please note; *It is the responsibility of the user to apply micro-droplets using thermoplastics or thermosets on the fibre when using the IFSS module.*

Delivering Measurement Solutions

e: sales@diastron.com

w: www.diastron.com

UK office 9 Focus Way | Andover | Hampshire | SP10 5NY | UK
t: +44 (0)1264 334700 | f: +44 (0)1264 334686

888 Sussex Boulevard | Broomall | PA19008 | USA
t: (610) 328-9038 | f: (610) 328-7157

US office

The company reserves the right to alter specifications and prices without giving prior notice.

© Dia-Stron Limited