

# Stand Out with Your Excellent Cosmetic Advice

## Multi Skin Test Center® MC 1000 & Software

moisture, oil, melanin, pH, elasticity & many more



Allows operating all C+K probes at the same time according to individual demands due to modular C+K probe technology.

C+K-measuring probes for following parameters can be connected: moisture, oil, pH, pigmentation & redness, transepidermal waterloss (TEWL – the barrier function of the skin), skin temperature and elasticity. Additionally the camera Visioscope® PC 35 can be connected.



The **Multi Skin Test Center® MC 750** offers the features like the MC 1000 except for the camera functions.

Please ask for our special transportation accessories.

**Technical Data:**  
 External power supply: Input: 100-240 V, 50-60 Hz, Output: 12 V, 4 A,  
 Dimensions: 46 x 22 x 9 cm; Weight without probes: approx. 1.8 kg,  
 Computer: Windows® 7 / Windows® 8 / Windows® 10; Interface: USB 2.0  
**Measuring principles:**  
 Moisture: capacitance (scale 0-99); Sebum: photometrical (scale 0-99); Elasticity: suction (scale 0-99); Melanin/Erythema: absorption/reflection (scale 0-99); pH: glass electrode (pH with one decimal); TEWL-Index: open chamber water flux measurement (scale 0-20)  
 Additional Visioscope® PC 35 or Visiopor® PP 34N

Technical changes may be made without prior notice.



[www.courage-khazaka.de](http://www.courage-khazaka.de)

Cologne, Germany

# Stand Out with Your Excellent Cosmetic Advice

## Multi Skin Test Center® MC 1000 & Software

moisture, oil, melanin, pH, elasticity & many more

### CSI-Software

- Very user-friendly and appealing software "Complete Skin Investigation" simplifying work with the system.
- Individual product recommendations for skin care, hair care and sun screen products which can be printed out for the customer.
- Easy increase of sales by offering a special service to distinguish from competitors and to secure customers retention.

### Software with 6 modules:

- Questionnaire combined with moisture, sebum and elasticity measurements as well as the camera measurements (optional) calculates **skin type** and recommends a set of skin care products (from cleansing to day and night care, additional care and food supplements) for women or men.
- Imaging** with Visioscope® PC35 and calculation of pores, spots and wrinkles. Easy comparison of customer images side by side.
- Sebum production & desquamation can be analyzed with the special foils Sebufix® and Corneofix®.

- Hair analysis** (measurement of sebum, colour, imaging for dandruff, hair thickness) with questionnaire and product recommendation.

- Answering questions and measurement of melanin content allows determination of most suitable **sun screen** and aftersun products, taking into account where, when and how long somebody is going to be exposed to the sun including the respective UV-index.

- Single measurements** (moisture, sebum, elasticity, pigmentation/erythema, pH, temperature and TEWL) in any order and on any body site to be saved and printed out.

- User can easily enter the products to be recommended in each programme according to the age of the person, the measurement results and answers to the questions.

- The appearance of the software can be easily changed by altering the name of the programme, exchange colours or the start screen picture, add a logo or even create a new questionnaire.

- Visioface® 1000 D** programme for analysis of full face photographs is integrated.

